# BAP

**EN - ES - FR**

**PRECISION SCALE CLASS II**

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V.1.2  
230421

## **EN**

<b>1. PRECAUTIONS BEFORE USING THE BALANCE</b>	<b>5</b>
<b>2. NAME AND TYPE OF THE INSTRUMENT</b>	<b>5</b>
<b>3. THE DISPLAY AND KEYPAD</b>	<b>5</b>
3.1 CONSTRUCTIONS	5
3.2 DISPLAY DESCRIPTION	5
3.3 KEYPAD FUNCTION	6
<b>4. FUNCTIONS AND OPERATION</b>	<b>7</b>
4.1 POWER-UP	7
4.2 WEIGHING MODE	7
4.2.1 DISPLAY RANGE	7
4.2.2 ZERO-SETTING	7
4.2.3 TARE FUNCTION	8
4.2.4 UNITS SELECTION	8
4.3 COUNTING MODE	8
4.4 PERCENT (%) WEIGHING MODE	9
4.5 DENSITY CALCULATION MODE (Optional)	10
<b>5. PARAMETERS SET-UP</b>	<b>13</b>
5.1 TO CHOOSE THE AUTO. OFF TIME	13
5.2 TO CHOOSE RS-232 TRANSMISSION METHOD	14
5.3 LABEL FORMAT	14
5.4 LABEL FORMAT	15
5.5 DATA BITS	15
5.6 TO CHOOSE U.WT RECOMPUTING	15
5.7 BACKLIGHT MODE	16
<b>6. RS – 232 INTERFACE</b>	<b>16</b>
6.1 MODE	16
6.2 FORMAT	16
6.3 DATA FORMAT	17
6.4 DATA FORMAT OF DENSITY CALCULATION	17
6.5 VARIABLES	17
<b>7. CALIBRATION PROCEDURE</b>	<b>18</b>
<b>8. ERROR MESSAGE</b>	<b>19</b>
<b>9. WEIGHT CONVERSION</b>	<b>19</b>

# ES

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## 1. PRECAUTIONS BEFORE USING THE BALANCE

For safe and dependable operation of this balance, please comply with the following safety precautions:

- Verify that the input voltage printed on the AC Adapter and the plug type matches the local AC power supply.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- Disconnect the scale from the power supply when cleaning the scale.
- Do not operate the scale in hazardous or unstable environments.
- Do not immerse the scale in water or other liquids.
- Do not drop loads on the platform.
- Use only approved accessories and peripherals, as available.
- Operate the scale only under ambient conditions specified in these instructions.
- Service should be performed by authorized personnel only.
- It must not be bumped against by other items or overloaded with excessively heavy weights (The load must not exceed the maximum capacity of the balance).
- Give a warm-up for 30mins before using the scale.

## 2. NAME AND TYPE OF THE INSTRUMENT

The BAP model weighing balance is a self-indicating weighing scale of Class II with single weighing range, an external AC mains adapter, and an internal rechargeable battery.

The balance consists of analogue to digital conversion, microprocessor control, power supply, keyboard, and a weight display contained within a single enclosure

## 3. THE DISPLAY AND KEYPAD


### 3.1 CONSTRUCTIONS

The balance is housed in a plastic enclosure. One display is placed on the front together with keyboard. It is intended to be used in shops as well as industrial environment.

A level indicator is fixed behind the enclosure and four adjustable leveling feet are located at the bottom

### 3.2 DISPLAY DESCRIPTION

The LCD display is described below

If the Capacity of the balance is 600.0g, the last digit will be the auxiliary display digit (a symbol indicating the digit):

If the Capacity of the balance is 300.00, 1200.0, 1500.0 and 3000.0g, the last digit will be blank in normal weighing mode only if the EXT.D key is pressed to extend the display.

Note: Once the balance is set to be in Trade mode, neither the auxiliary display nor the extended display function can work.

### 3.3 KEYPAD FUNCTION

KEY	DESCRIPTION	
	PRIMARY FUNCTION	SECONDARY FUNCTION
	To power on or power off the balance.	
	To choose weighing, parts counting, percentage weighing or density function mode.	
	To extend the display.	To take samples in parts counting mode or percentage weighing mode.
	To select the available units.	To increase value or page up
	To subtract the container's weight	To decrease value or page down

	To return the display to "0"	To return to last step in setting
	To select the backlight type	To move the digit to right
	To transmit or print the data	To confirm the current setting in setting mode

## 4. FUNCTIONS AND OPERATION

The balance can be available for operation from mains at 110~240V (50/60Hz) with an external AC/DC adapter with 12V output voltage and an optional internal 6V rechargeable battery.

### 4.1 POWER-UP

When power up the balance, the display will first show the software version for 2 seconds and then perform a self display test. Then it will show the counting times for about 1 second. After that it will automatically establish the current weight as a new zero reference point.

### 4.2 WEIGHING MODE

#### 4.2.1 DISPLAY RANGE

The display can show weight Tare weight to Max+9e (Gross weight).

#### 4.2.2 ZERO-SETTING

Zero-setting range:  $\pm 2\%$  of Max

Zero-setting is only available when the load receptor is not in motion (The display is stable).

#### SEMI-AUTOMATIC ZERO-SETTING

When there is a minor weight displayed (without load on the pan), press ZERO key to return to the display to zero, then the Zero symbol appears which indicates that the display now is at the zero point.

#### ZERO-TRACKING

The instrument is equipped with a zero-tracking feature which operates at a speed of  $\pm 0.25e/s$  and only when the indicator is at gross zero and there is no motion in the weight display.

#### INITIAL ZERO-SETTING


Initial zero-setting range:  $\pm 10\%$  of Max

When the initial zero-setting range is below  $-10\%$  of Max, the display shows "-----" and the buzzer sounds. To return the display to normal mode, please adding more weights until the initial zero-setting range is between  $\pm 10\%$  of Max.


When the initial zero-setting range is beyond  $+10\%$  of Max, the display shows "-----" and the buzzer sounds. To return the display to normal mode, please take off weights until the initial zero-setting range is between  $\pm 10\%$  of Max.

### 4.2.3 TARE FUNCTION

Press TARE key to subtract the current displayed weight value as the tare weight value (Consecutive tare operations are permitted). The NET symbol appears and the display will show a net weight of the object to be weighted.

Then take off the container, the display will show a minus "-" weight value which equals to the container's weight.

### ELIMINATE THE TARE WEIGHT

To eliminate the tare weight by pressing  Key with no object on the scale.

Note:


1. When performing the tare operation(s), make sure the instrument is in stable.
2. If consecutive tare operations have been performed, to clear the tare value, please remove all the containers on the pan first and then press  key.
3. Tare range is the full capacity of balance for all models.

### 4.2.4 UNITS SELECTION

Press  key to choose the weighing units and the display will be changed to the new value with the units shown. There are up to 13 units of weight that can be enabled. Refer to Parameter Set-up Section.

### 4.3 COUNTING MODE

In normal weighing mode, press  key to enter into counting mode, and the display shows the sample size:


Press the  key to select the sample size from 10, 20, 50, 100, 200, 500 y 1000.

### SAMPLING

Put on some samples with a sample size the same as the current displaying one. Then press the  key and the displayed value starts flickering. Several seconds later, the display shows as below:

Take off the samples on the pan. The sampling operation is finished. Then the balance can be used for counting.


Press  key to return to normal weighing mode (Pressing  key after it returns to weighing mode will bring back to counting mode and the unit weight data is still stored for use).

**Note:**

1. If a container needs to be used, press TARE key to subtract its weight before entering into the counting mode.
2. If the unit weight is too small (less than 0.2 of a scale division), the display will show: --or-- PCS.

## 4.4 PERCENT (%) WEIGHING MODE

In normal weighing mode, press the  key twice to enter into percent weighing mode.


Press  or  key to change to "100.00%".

### SAMPLING

Put the item to be considered 100% on the pan and press the  key. Then the displayed value starts flickering. Several seconds later, it stops flickering.

Take off all the weight on the pan and the display shows as below:

Then the sampling operation is finished. The balance can be used for percent weighing.

Press  key to return to normal weighing mode. (Pressing  key after it returns to weighing mode will bring back to percent weighing mode and the unit weight data is still stored for use).

**Note:**

1. If a container needs to be used, press  key to subtract its weight before entering into the counting mode.
2. If the unit weight is too small (less than 0.2 of a scale division), the display will show: --or-- PCS.
3. The memory of sampling size data will be cleared automatically every time the  key is used to select parts counting or percent weighing.

#### 4.5 DENSITY CALCULATION MODE (Optional)


The balance can calculate sample's density through measuring sample weight in air, weight in water and water's density in current temperature

**PROCEDURE**


Press the  key 3 times in normal weighing mode to switch to density calculation mode. Then the display shows a temperature value.

The balance can be equipped with an optional hanging hook which can be used in Density Calculation mode.

1. Open the rubber seal on the bottom of the balance.

2. Install the hanging hook.

**A. CHOOSE WATER'S TEMPERATURE.**


Use the thermometer to test water's temperature in current situation. Use

 or  key to choose the corresponding temperature value. Press  ENTER key to confirm.

B. Then the display shows the water's density under the temperature you have chosen. And press  key to confirm.

For example: 0.99705 $\rho$ -L  
 $\rho$ -L: Indicate water density

**C. SAMPLE WEIGHT IN AIR**

Then the display shows "0.00g in-A" Hang the sample on the hook (If a container is used, press  key to subtract the container's weight first), the display will show the sample's weight in air.

In-A: Indicate weight in air  
 The display shows as below


Press  key to confirm

**D. SAMPLE WEIGHT IN WATER**

Then the display shows “0.00g in-L” .Hang the sample on the hook (If a container is used, hang the container on the hook and put it into the water, then press TARE key to subtract the container’s weight first), and then put the sample in the water. The display will show a weight value, which indicates the sample’s weight in water.

in-L: Indicate weight in water

The display shows as below

Press  key to confirm

**E.** Then the display will show a value which indicates the density of the sample. The density calculation is finished. ( ρ-O: density of the sample)

For example: “33.6841ρ-O”

$$\rho-O = \frac{100.00}{100.00-97.04} \times 0.99705$$

Press  key to confirm and do the density calculation of next object.

Press  key to exit from the density calculation mode and return to weighing mode.

**REMARKS**

1. The density unit is g/cm<sup>3</sup>, so if the weighing unit is not “g” when entering into this mode, it will automatically be changed to “g.” And when exiting from this mode, the weighing unit will be back to the original one.
2. Press  key to exit from this mode at any time.
3. Zeroing, Tarring and cancelling tare functions are available in this mode.
4. Density calculation formula:

$$\rho-O = \frac{W_A}{(W_A - W_L)} \times \rho-L$$

ρ-O	Density of the sample
ρ-L	Density of the water
W <sub>A</sub>	Sample weight in air
W <sub>L</sub>	Sample weight in water

5. Press  key in step “e” to do the density calculation for next object, the water temperature will be the same as the last one. But if exiting from this mode, the water temperature will be the default one when entering into the mode for the next time

## 6. Water's density in different temperature

Tem. (°C)	Den. (g/cm <sup>3</sup> )	Tem. (°C)	Den. (g/cm <sup>3</sup> )	Tem. (°C)	Den. (g/cm <sup>3</sup> )	Tem. (°C)	Den. (g/cm <sup>3</sup> )
0	0.99984						
1	0.99990	11	0.99961	21	0.99799	31	0.99534
2	0.99994	12	0.99950	22	0.99777	32	0.99503
3	0.99996	13	0.99938	23	0.99754	33	0.99471
4	0.99997	14	0.99925	24	0.99730	34	0.99438
5	0.99996	15	0.99910	25	0.99705	35	0.99404
6	0.99994	16	0.99894	26	0.99679	36	0.99369
7	0.99990	17	0.99878	27	0.99652	37	0.99333
8	0.99985	18	0.99860	28	0.99624	38	0.99297
9	0.99978	19	0.99841	29	0.99595	39	0.99260
10	0.99970	20	0.99821	30	0.99565	40	0.99222


## 5. PARAMETERS SET-UP

The parameter set-up menus allow the balance to be set as required by the user.

Press the  and  key at the same time, and hold for 2 seconds, then the display will enter into parameter setting mode.

### 5.1 TO CHOOSE THE AUTO. OFF TIME

The display show as below:


Press  or  key to move to last setting or next setting; or press  key to enter into current setting.

Press the  or  key to choose the Auto. off time from 2m, 5m, 8m or no (Disabled to auto. off). (Default: no)

Then press the  key to confirm and move to next setting or press the  key to move to next setting.

## 5.2 TO CHOOSE RS-232 TRANSMISSION METHOD

The display show as below:

Press  or  key to move to last setting or next setting; or press  key to enter into current setting.

Press the  or  key to choose transmission method.

(Default: Close)

tr.m.LP = used for LP50 printer (pressing  key)

Tr.A.PC

Tr.A.dT

Tr.A.LP

Tr.disp =used for connection with extra display

Tr.m.PC = used for PC ((pressing  key)

Tr.m.dt = used for GODEX BP DT-4 printer ((pressing  key)

Then press the  key to confirm and move to next setting or press the  key to move to next setting.

## 5.3 LABEL FORMAT

The display show as below:

Press  or  key to move to last setting or next setting; or press  key to enter into current setting.

Use BL key to move the digit; use the  or  key to key in the digits: FORM00-FORM99. (Default:

FORM00)

Then press the  key to confirm and move to next setting or press the  key to move to next setting.

#### 5.4 LABEL FORMAT

The display show as below:

# rate

Press  or  key to move to last setting or next setting; or press  key to enter into current setting.

Press the  or  key to choose the desired baud rate: “1200”, “2400”, “4800” or “9600”. (Default: 9600)

Then press the  key to confirm and move to next setting or press the  key to move to next setting.

#### 5.5 DATA BITS

The display show as below:

# Data b

Press  or  key to move to last setting or next setting; or press  key to enter into current setting.

Press the  or  key to choose the data bits to be 7

(bit.e7) or 8 (bit.n8)). (Default: 8)

bit.e7: data bits 7, even parity bit

bit.n8: data bits 8, no parity bit

Then press the  key to confirm and move to next setting or press the  key to move to next setting.

#### 5.6 TO CHOOSE U.WT RECOMPUTING

The display show as below:

# AC Y

Press  or  key to move to last setting or next setting; or press  key to enter into current setting.

Press the  or  key to choose to be enabling (y) or disabling (n). (Default: n)

Then press the  key to confirm and move to next setting or press the  key to move to next setting

## 5.7 BACKLIGHT MODE

The display show as AUTO/ OFF / ON

Press  or  key to move to last setting or next setting; or press  key to enter into current setting.

Press the  or  key to choose the backlight mode: AUTO, OFF and ON. (Default: AUTO)

AUTO: The backlight will be on automatically when a load exceeds 10e on the pan or pressing a key.

ON---The backlight will always be on.

### Note:

1. When the scale is in the parameter setting mode, pressing  key means the selected setting value is confirmed.

2. When the scale is in the parameter setting mode (any type of setting mode), pressing  key will return to normal weighing mode without making a change.

## 6. RS – 232 INTERFACE

### 6.1 MODE

Standard EIA-RS-232 C's UART signal

### 6.2 FORMAT

Standard EIA-RS-232 C's UART signal

*Baud rate:* 1200, 2400, 4800, 9600

*Data bits* 7/8 BITS

*Stop bit* 1 BIT

*Code* ASCII

RS-232 connector is a 9 pin D-subminiature socket.

*Input Pin* 2

Output Pin 3  
Signal Ground Pin 5

### 6.3 DATA FORMAT

HEAD 1, HEAD 2 , DATA UNIT <CR><LF>  
1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20

HEAD1 (2 BYTES)

HEAD2 (2 BYTES)

OL - Overload, under load  
ST - Display is stable

NT – NET Mode  
US - Display is unstable

DATA (8 BYTE)

2C (HEX) = “,” (COMMA)                    0D (HEX) = <CR>  
2D (HEX) = “-” (MINUS)                0A (HEX) = <LF>  
2E (HEX) = “.” (DECIMAL POINT)  
20 (HEX) = “ ” (SPACE)  
UNIT (4 BYTE), text for units selected, may include leading spaces.

#### EXAMPLES:

EX: +0.876g, when it is stable and net value as:  
ST, NT, + 0.876 g <CR><LF>  
EX: -1,568lb, when it is unstable and net value as:  
US, NT, - 1.586 lb <CR><LF>

### 6.4 DATA FORMAT OF DENSITY CALCULATION

It is available when the transmission method is set to be “Tr.PC”

<p><math>\rho-L: 0.99705 \text{ g/cm}^3</math></p>	Press  key to confirm waterdensity.
<p><math>W_A: 127.19\text{g}</math></p>	Press  key to confirm sample weight in air.
<p><math>W_L: 92.87 \text{ g}</math></p> <p>.....</p>	Press  key to confirm sample weight in water.
<p><math>\rho-O: 3.69615 \text{ g/cm}^3</math></p>	Press  key to confirm sample weight in density.

$\rho-O$ -----Density of the sample     $\rho-O$ -----Density of the water  
 $W_A$ : -----Sample weight in air     $W_L$ -----Sample weight in water


### 6.5 VARIABLES

NO.	VARIABLE NAME		SPECIFICATION	APPLICATIONMODE	SIZE
	LP-50	DT			
1	NWA	V0	Net weight (with “. ”)	W/C/%	7
2	NWB	V1	Net weight	W/C/%	6
3	TWA	V2	Tare weight (with “. ”)	W/C/%	7

4	TWB	V3	Tare weight	W/C/%	6
5	GWA	V4	Gross weight (with ".")	W/C/%	7
6	GWB	V5	Gross weight	W/C/%	6
7	QUA	V6	Quantity/Percentage (.)	C/%	7
8	QUB	V7	Quantity/Percentage	C/%	6
9	UWA	V8	Unit weight (with ".")	C	7
10	UWB	V9	Unit weight	C	6
11	UNT	V10	Weighing Unit	W/C/%	3

## 7. CALIBRATION PROCEDURE

Switch on the balance, and press the CAL switch during the self-test.

When the self-test procedure finished, the display will show as below:

CAL

When it shows "CAL," press  key to start calibration procedure.

0.00 g

- Touch the pan lightly by hand, and several seconds later, once it take stability again, the display shows the weight value of first calibration mass (1/3 of full capacity). Place the displayed weigh on the pan.

100.00 g

- Several seconds later, once it take stability, the display shows the weight value of second calibration mass (2/3 of full capacity) Place the displayed weigh on the pan.

<h1 style="margin: 0;">200.00 g</h1>
--------------------------------------

- Several seconds later, once it take stability, the display shows the weight value of third calibration mass (full capacity) Place the displayed weigh on the pan.

<h1 style="margin: 0;">300.00 g</h1>
--------------------------------------

Several seconds later, the buzzer tweets. Then the calibration procedure is finished. Switch off the balance & take away all the weights; then switch on again and the balance can be used.

**Note:** -If during the loading of calibration mass process scale arrive to take stability whithout have the properly mass on the pan, it will assume stable weigh as asked weigh, producing a wrong calibration. You must press ESC key for to start calibration process again.

When balance is on, the display may show E1 or E2, which can be eliminated by recalibrating the balance.

CODES	CAUSES
E1	The Span calibration data loses
E2	The initial zero drifts when power on.

Before calibrating the balance, it must have been turned on for at least 30 min in the room where the calibration is to be carried out.

## 8. ERROR MESSAGE

During operation, the scale may show some error messages.

CODES	POSSIBLE CAUSES
E1	EPROM data loss or MCU is damaged
E2	The zero of the load cell is out of the initial zero range.
E5	The sample weight is not enough
OL	Overload

Correct your operation and contact the dealer or manufactory for further assistance.

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